

PRODUCT SPECIFICATION SHEET


HAP STUD MOUNTING BRACE

Small HAP Stud Mounting Brace

Model Number: PL-1111-MB1

Large HAP Stud Mounting Brace

Model Number: PL-1111-MB10




Model PL-1111-MB1

The Small Mounting Brace can be used with	All Small, Medium and Large HAP Hangers
Small Brace Adjusts from	10-3/8" to 18"
Material	Galvanized Steel
Max Load	DNA
Length	10-3/8" to 18"
Width	1-1/4"
Height	5/8"
Metal Detection	Yes
Finish	Galvanized
Inside or Outside wall Installation	Determined by required pipe location

Model PL-1111-MB10

The Large Mounting Brace can be used with	All Small, Medium and Large HAP Hangers
Large Brace Adjusts from	15" to 26"
Material	Galvanized Steel
Max Load	DNA
Length	15" to 26"
Width	1-1/4"
Height	5/8"
Metal Detection	Yes
Finish	Galvanized
Inside or Outside wall Installation	Determined by required pipe location

- The brace length is adjustable to fit your specific needs.
- Locator dimples allow exact positioning of a single or multiple HAP Hangers.
- Dimples allows screws to start quickly and prevents stripping out when attaching HAP Hangers.
- The mounting legs are bendable and can be used on metal or wood studs.



Product Submittal	
Date:	
Job Name:	
Model Number	Quantity
Architect / Owner:	
Contractor:	
Notes	